


# Foley Catheter



Code: UR-07

Latex, Siliconised

Smooth tapered tip of catheter facilitates easy introducing into the urethra

The drainage eyes are accurately formed to permit effective drainage

The symmetrical balloon expands equally in all directions so as to perform its function of retaining bladder efficiently

The smooth outer surface is specially treated with medical grade silicone fluid which facilitates easy passage through urethra

Colour Coded sleeves are useful for easy and rapid size identification

Available with Plastic or Rubber valve

Available Sizes: 6 to 24

## Options Available

- \* 100% silicon
- \* 2 Way / 3 Way
- \* Child / Adult Size
- \* Poly Pack or Paper-Poly Pack