

## Suction Catheter

Code: AS-01

Designed for tracheal / bronchial suction. Smooth clear PVC ensures efficient mucus suction while minimizing potential damage to the delicate tracheal mucous membranes.

Specially designed for easy passage through plastic endotracheal and tracheotomy tubes with minimum irritation to patient.

Open distal end with elliptical eye for easy and efficient suction.


Color-coded connector at proximal end for instant identification of catheter size.

Provided with X-ray opaque line for radiological identification.

Size (FG) : 8, 10, 12, 14, 16, 18, 20, 22, 24

### Options Available

- \* Various sizes and lengths
- \* With or without XRO Line
- \* Poly Pack or Paper-Poly Pack
- \* Inner Bag or Inner Box



## Suction Catheter with Thumb Controller

Code: AS-02

Manufactured from non toxic non irritant medical grade PVC

Suction catheters are suitable for removal of secretion from mouth oropharynx trachea and bronchial tubes

Distal eye is open with one lateral eye

Provided with " T " type vacuum control valve with thumb control for proper maneuvering

Frozen surface tubing for smooth intubation

Colour Coded for instant size identification