

Nebulizer Mask with Tubing



Code: AS-11

Options Available

- * ADULT / PAEDIATRIC
- * Poly Pack or Paper-Poly Pack
- * Inner Bag or Inner Box

Made from special transparent non toxic material.

Nebuliser chamber capacity 6/7cc.

Facemask with aluminium clip designed to fit comfortably over the patient's mouth.

The kit contains Jet Nebuliser, oxygen mask with nose clip, elastic band and 180 cm long tubing .

The connecting tube has internal serrations (STAR LUMEN).

ETO Sterile and Pyrogen free