

## Nasal Oxygen Catheter



Code: AS-03

Designed for administration of oxygen.

Green colour catheter made of implantation tested PVC having open distal end with multiple eyes to facilitate dispersion of oxygen and prevent oxygen burn.

Colour coded connector for size identification.

Sterile, non-toxic & pyrogen free.

Available with 8 holes

Size: 6, 8, 10, 12, 14.

### Options Available

- \* Various sizes and lengths
- \* With or without XRO Line
- \* Poly Pack or Paper-Poly Pack
- \* Inner Bag or Inner Box