


Tracheostomy Tube Plain

Code: AS-10

Designed for long-term tracheostomy
Manufactured from non-toxic siliconised
PVC to protect delicate mucosal tissues
Provided with soft flexible flange for
easy fixation
Terminal end is suitable for 15mm
connector
Individual Packing in Blister Pack
ETO Sterile, Latex free
Available sizes: 4.0mm to 10mm
Paper Poly Packing



Tracheostomy Tube Cuffed

Code: AS-09

Non-return valve with indicator
facilitate the cuff inflation and
deflation

Thin, smooth cuff inflates
symmetrically eliminating any risk of
herniation

Terminal end is suitable for 15mm
connector

Individual Packing in Blister Pack

ETO Sterile, Latex free

Available sizes: 5.0mm to 10mm