



Code: AS-07A & 08A

Re-Inforced Endotracheal Tube (Plain & Cuffed)

Re-Inforced Endotracheal Tubes permit complete airway control for short term intubations.

The tubes are firm enough for ease of intubations but soften at body temperature to conform to the anatomy of the upper respiratory track.

Kink-resistant, minimize possibility of tube occlusion.

Oral/Nasal version a flexible approach to the route of intubations.

One-way valve for efficient and easy inflation and deflation of cuff.

Comprehensive size range to suit all adolescent and adult tracheal sizes.

EO Sterile, Latex Free

Individual Packing in Blister Packing

Available Sizes:

With Cuff: 3.0 to 10 sizes

Without Cuff (Plain): 2.0 to 10 sizes

Size of the tube clearly marked on the tube

Cuff is optimally designed for routine intubation

One way valve for efficient and easy inflation and deflation of cuff

Murphy eye, high volume, low pressure cuff, oral/nasal

ETO Sterile, Latex free