


# Umbilical Catheter



Code: MS-05

Designed for intermittent or continuous access to the umbilical artery or vein of newly born or premature infant.

Open distal end without lateral eyes to eliminate chance of clot forming in the blind spaces.

Proximal end fitted with luer mount for convenient connection to I.V. Set.

Sterile, non-toxic and pyrogen free.

Size (FG): 3, 4, 5, 6, 8, 10.

## Options Available

- \* Various sizes and lengths
- \* With or without XRO Line
- \* Poly Pack or Paper-Poly Pack
- \* Inner Bag or Inner Box