

Umbilical Cord Clamp



Code: MS-06 & 07

Designed for clamping the umbilical cord soon after the birth.

Provided with double purpose security lock with a click to indicate correct locking and protect against accidental re-opening.

Provided with grooves all along the length to prevent slip of the umbilical cord and to retain it in same position.

Clamp arms inclined at a suitable angle for easy movement for clamping.

Finger grip to ensure safe and convenient handling, particularly for wet gloves.

Available with re-opening facilities.

Sterile, non-toxic and pyrogen free.

Options Available

- * Poly Pack or Paper-Poly Pack
- * Inner Bag or Inner Box