


Silicon Foley Catheter



Code: UR-07A

Made 100% medical grade silicon, have good biological compatibility

Smooth tapered tip of catheter facilitates easy introducing into the urethra

The drainage eyes are accurately formed to permit effective drainage

The symmetrical balloon expands equally in all directions so as to perform its function of retaining bladder efficiently

The smooth outer surface is specially treated with medical grade silicone fluid which facilitates easy passage through urethra

Colour Coded sleeves are useful for easy and rapid size identification

Different check valve include the type with/without Luer Lock

Available Sizes: 8 to 24

Options Available

- * Latex with Silicon coated
- * 2 Way / 3 Way
- * Child / Adult Sizes
- * Poly Pack or Paper-Poly Pack