

# Enteral Feeding Bag



Code: GFB-01 & PFB-01

## Gravity / Pump Systems

- 1000ml capacity graduated bag with long tubing
- Rigid neck for easy filling and handling
- Large Easy Cap Opening
- Strong, dependable hanging ring
- Easy to read graduations and easy view translucent bag material
- Bottom exit port allows complete drainage
- Non-I.V. compatible connector
- Sterile
- Single packed



Code: EFT-01

**Enteral Feeding Set also available**

## Important Safety Information

- \* Assure that set is connected only to enteral feeding tube.
- \* To avoid product contamination problems, this set should be replaced at least every 24hrs., or as needed.
- \* Single patient use only
- \* See individual product labeling for Directions for Use and important Safety Information