



Code: AS-20

SPIROMETER

Spirometer is lung exerciser. Most commonly used after surgery. It encourages deep and prolonged inspiration and helps patients in faster recovery post operatively.

Specifications:

- ◆ Three color-coded balls in each chamber rise proportionately to the inspiration exerted by the user provide a visual incentive for the patient.
- ◆ Connecting tube with 12mm OD connector and mouth piece.
- ◆ Flow rates 600ml/sec, 900ml/sec, and 1200ml/sec by using different colors of ball for easy to identify the flow rate.
- ◆ Dis-assembled into parts for cleaning and disinfection.
- ◆ Non sterile, individually packed in a box.
- ◆ Material: ABS for chamber, PE for ball and mouthpiece, and EVA for tubing.

