

Anesthesia Face Mask



Code: AS-26

- ◆ Transparent cushioned, crystal clear, flexible 22 F ISO connection.
- ◆ Air filled cushion to provide safety and comfortable use.
- ◆ Suitable for use in respiratory care of anesthesia and resuscitator
- ◆ Body with anti-static property & good grip.
- ◆ Inflatable cushion and rounded edges offer optimal comfort.
- ◆ Removable color-coded head strap hook rings. Allows use of traditional head straps. Removable for use in hand-held procedures.
- ◆ Available Sizes : 0 - Neonatal , 1 - Infant , 2 - Pediatric , 3 - Adult Small , 4 - Adult Medium , 5 - Adult Large.
- ◆ Individually packed in Poly Pack.