




## BURETTE SETS

- Pediatric Burettes (Volumetric): Available in 100ml and 150ml.
- 60 Drops / ml.
- Burette chamber is scaled for accurate measurement of medication for critical and pediatric /neonatal care.
- Options available: Floating Valve, Fluid Filter, Y-Injection Site, Flash Bulb etc.
- ISO 8536-4 compliant spike ISO 594/1 & ISO 594/2 compliant end-connection
- Available with Latex and Latex Free.
- Soft Clear Flexible Tubing
- Sterile, Non-Toxic and Pyrogen-Free.
- Sterile Pack - ETO (Ethylene Oxide) Sterilization validated according to Harmonized Standards


### Available Specifications



**Code: SP-05**

81" Burette Set 150ml


Vented Spike 60 Drops-ml with Handle/Latex Free MV Cylinder/ Shut-off Valve/ Chamber/Fluid Filter - 15micron/Slide Clamp/Flow Regulator/Latex Free 'Y' Site/Rotating Male Luer Lock Adaptor with Cap/2.7mm x 3.9mm Tubing-Length 25+135+15cm/ EO Sterile/Paper Poly Pack/ Quantities per Carton = 80 pcs



**Code: IT-27-A**

72"/70" Burette Set 150/100ml


Non-Vented Spike 60 Drops-ml with Handle/MV Cylinder with Shut-Off Valve/ Chamber/Fluid Filter/Slide Clamp/Flow Regulator/Bulb Latex/Luer Tip Adaptor/ Needle 23g x 1" with Protector/2.7mm x 3.9mm Tubing-Length 25+130cm/ EO Sterile/Poly Pack/ Quantities per Carton = 120 pcs



**Code: IT-27-B**

73"/71" Burette Set 150/100ml

Non-Vented Spike 60 Drops-ml with Handle/MV Cylinder with Shut-Off Valve/ Blood Chamber/ Blood Filter/80mm Drip Chamber/Slide Clamp/Flow Regulator/ Bulb Latex/Luer Tip Adaptor/Needle 23g x 1" with Protector/2.7mm x 3.9mm Tubing Length 25+130cm/EO Sterile/Poly Pack/ Quantities per Carton = 120 pcs



**Code: IT-27-C**

73"/71" Burette Set 100ml/150ml

Vented Spike 60 Drops-ml with Handle/MV Cylinder with Shut off valve/ Fluid Filter/Drip Chamber/Big Flow Regulator/Latex Free 'Y' injection site/ Luer Lock with Protective Cap/PVC Tube 2.7x3.9mm, Length 25+130cm, EO Sterile/Poly Pouch Quantities per Carton = 120 pcs