

C.V.P. MANOMETER



Code: IT-37

Infusion Therapy

Suitable for continuous or intermittent monitoring of central venous pressure during infusion.

Manometer tube with graduation from -4cm to +34cm and attached to three way connector.

125 cm long extension tube is fitted with female luer lock at one end and male luer lock on the other end for safe connection to C.V.P. Catheter.

Latex rubber flash ball is provided for extra medication.

Sliding indicator provided to record previous reading.

Moulded clamps with latex rubber strings are provided for easy fixing on I.V. Stand.

Individual packing.

ETO Sterile and Pyrogen free

CVP Manometer