

## IV Set with Precision Flow Controller (Dial type)



Ideally suited for gravity infusions  
 Available in both vented (for bottles and bags) and non-vented (for Plastic Bottle bottles) sharp spike for better penetration.  
 ISO 8536-4 compliant spike ISO 594/1 & ISO 594/2 compliant end-connection  
 Precision Flow Controller (Dial type) can be set from 5-250ml/hour.  
 Latex Free  
 Soft Clear Flexible Tubing  
 Sterile, Non-Toxic and Pyrogen-Free.  
 Sterilized by E.O. (Ethylene Oxide) Gas.

Infusion Sets - Dial type

### Available Specifications



**Code: IT-20-A**

83" Vented Spike with Guard/20 drops-ml/Drip Chamber 60mm/Fluid Filter - 15micron/Flow Regulator/  
 Dial Flow Regulator/ Latex Free 'Y' Site /Rotating Luer Lock Male Adaptor/  
 2.7mm x3.9mm PVC Tubing-Length 50+130+15cm/EO Sterile/Paper Poly Pack  
 Quantities per Carton = 240 pcs



**Code: IT-20-B**

83" Vented Spike with Guard/60 drops-ml/Drip Chamber 60mm/Fluid Filter - 15micron/Micro Drip  
 Tube/Flow Regulator/Dial Flow Regulator/ Latex Free 'Y' Site /Rotating Luer Lock Male  
 Adaptor/2.7mm x3.9mm PVC Tubing-Length 50+130+15cm/EO Sterile/Paper Poly Pack  
 Quantities per Carton = 240 pcs