

Oxygen Mask



Code: AS-05A

- ◆ Designed for convenient oxygen therapy.
- ◆ Molded face mask with adjustable elastic strap for proper position of mask on the mouth and nasal area.
- ◆ Provide with 2-meter long star lumen tube to ensure continuous flow of oxygen.
- ◆ Proximal end fitted with soft funnel shape connector for easy connection to oxygen source.
- ◆ Sterile, non-toxic and pyrogen-free .

NON-REBREATHING MASK

Code: AS-05B



Options Available

- * Adult/Paediatric/Infant Sizes
- * Poly Pack or Paper-Poly Pack
- * Inner Bag or Inner Box