


# BT Sets - Economy Models

Transfusion Therapy - Economy Models




- Sharp spike vented (for bottles & bags) and non-vented (for plastic bottles) for better penetration.
- ISO 8536-4 compliant spike, ISO 594/1 & ISO 594/2 compliant end-connection
- Superior quality of latex tube for self-sealing.
- Specially designed roller clamp for better flow control.
- Top quality needle.
- Available with or without air-vent
- Non-Toxic, Pyrogen Free & Sterile by EO

### Available Specifications



**Code: TT-01**

63" Bottle Needle Blood Administration Set Non Vented  
 16g x 1.25" Bottle Needle Single Slot with Cap/20 drops-ml/Drip Chamber 80mm/Blood Filter/  
 Tube Latex/Flow Regulator/Luer Tip Male Adaptor /Needle 18g x 1.5" with Protector/  
 2.8mm x 3.8mm PVC Tubing- Length 135cm/Extra Airway Needle 16g x 1.25" Single Slot  
 Needle with Cotton and 5cm long Tube/EO Sterile/Poly Pack  
 Quantities per Carton = 500 pcs



**Code: TT-02**

60" Blood Administration Set Non Vented  
 Non-Vented 20 drops-ml Spike/Drip Chamber 80mm/Blood Filter/Tube Latex/Flow Regulator/  
 Luer Tip Male Adaptor /Needle 18g x 1.25" with Protector/2.8mm x 3.8mm  
 PVC Tubing- Length 135cm/Extra Air Spike with Protector/EO Sterile/Poly Pack  
 Quantities per Carton = 500 pcs