



Extension Lines

- Made from non-toxic medical grade PVC compound
- Individually pressure tested
- Available in assorted length
- 100% Leak proof and no blockage
- ISO 594/1 & ISO 594/2 compliant end-connection
- Options available with Precision Flow Controller or 3 Way Stopcock
- Soft Clear Flexible Tubing
- Sterile, Non-Toxic and Pyrogen free.
- Sterilized by E.O. (Ethylene Oxide) Gas.

Available Specifications



Code: IT-34

Extension Set with Precision Flow Controller

1 No. Male and 1 No. Female Luer Lock Adaptor with Cap/Latex Free 'Y' Site/Precision Flow Controller 5 to 250ml-hr/ 2.8 x 4.1 PVC Tubing-Length 100cm/Sterile/Paper Poly Pack
Quantities per Carton = 500 pcs



Code: IT-35A

Extension Set High Pressure/Pressure Monitoring

1 No. Male and 1 No. Female Luer Lock Adaptor with Cap/1.4 x 2.8 PVC Tubing/ Sterile/Paper Poly Pack/ Length (cm) - 10, 25, 50, 100, 150 & 200
Priming Volume (ml) - 0.16, 0.39, 0.78, 1.60, 2.40 & 3.10 (approx.)
Quantities per Carton = 600 pcs



Code: IT-35B

Extension Set with High Capacity

1 No. Male and 1 No. Female Luer Lock Adaptor with Cap/2.8 x 4.1 PVC Tubing/ Sterile/Paper Poly Pack/ Length (cm) - 10, 25, 50, 75, 100, 150 & 200
Priming Volume (ml) - 0.65, 1.55, 3.10, 4.68, 6.20, 9.30 & 12.40 (approx.)
Quantities per Carton = 600 pcs



Code: IT-36A

Extension Set with 3 Way Stopcock


1 no. Male Luer Lock Adaptor with Cap/1 no. 3 Way Stopcock with 2 nos. Female Luer Lock Adaptors with Cap/2.8 x 4.1 PVC Tubing/Sterile/Paper Poly Pack/
Length (cm) - 10, 25, 50, 75, 100, 125, 150 & 200
Priming Volume (ml) - 0.65, 1.55, 3.10, 4.68, 6.20, 7.75, 9.30 & 12.40 (approx.)
Quantities per Carton = 500 pcs



Code: IT-36B

Extension Set - Mini Volume with T Connector

T Connector with Injection Site and Male Luer Lock Adaptor with Cap/Slide Clamp/ Female Luer Lock Adaptor with Cap/1.4 x 2.8mm PVC Tubing -Length 18cm/ Sterile/Paper Poly Pack/ Quantities per Carton = 500 pcs



Code: IT-36C

Bifurcate/Trifurcate Extension Set

1 no. Male Luer Lock Adaptor with Cap/ 2 or 3 nos. Female Luer Lock Adaptors with Caps/Connectors/2 nos. Slide Clamp/ 2.8 x 4.1mm PVC Tubing -Length 10cm + 2x25cm/ Sterile/Paper Poly Pack/ Quantities per Carton = 700 pcs

