



Code: GE-01

Feeding Tube

- ◆ Designed for tiny infant and neonate for feeding.
- ◆ Made from implantation-tested PVC, non irritant to delicate mucous membranes
- ◆ Super smooth coned tip having two lateral eyes at distal end.
- ◆ Color coded female luer mount provided for safe feeding with syringe or funnel.
- ◆ A radio opaque strip co-extruded as an integral part of the catheters wall.
- ◆ Sterile, non - toxic and pyrogen free.
- ◆ Size (FG) : 3, 4, 5, 6, 7, 7.5, 8, 9, 10, 11, 12

Options Available

- * Various sizes and lengths
- * With or without XRO Line
- * Poly Pack or Blister Pack
- * Inner Bag or Inner Box