



## Infusion Sets for Infusion Pump

Infusion Therapy


- Available in vented (for bottles and bags) and non-vented (for Plastic Bottle bottles) Sharp spike for better penetration
- ISO 8536-4 compliant spike ISO 594/1 & ISO 594/2 compliant end-connection
- Options available for Latex Free
- Soft Clear Flexible Tubing
- Sterile, Non-Toxic and Pyrogen-Free.
- Sterile Pack - ETO (Ethylene Oxide) Sterilization validate according to Harmonized Standards

### Available Specifications



Code: SP-03-A

62" IV Administration Set Non-Vented  
 Non-Vented Spike 10 Drops-ml with Guard/Drip Chamber 60mm/Fluid Filter-15 micron/  
 Latex Free 'Y' Site/Flow Regulator/Luer Tip Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-  
 Length 135+15cm/EO Sterile/Paper Poly Pack/  
 Quantities per Carton = 300 pcs



Code: SP-03-B

62" IV Administration Set Vented  
 Vented Spike 20 Drops-ml with Guard/Drip Chamber 70mm/Fluid Filter-15 micron/  
 Latex Free 'Y' Site/Flow Regulator/Luer Lock Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-  
 Length 135+15cm/EO Sterile/Paper Poly Pack/  
 Quantities per Carton = 300 pcs

For Infusion Pump