

## Heat Moisture Exchange Filter (HME Filter)



Code : AS-22

- ◆ The Heat and Exchange Moisturizers Filter (HMEF) is used in mechanically ventilated patients intended to help for preventing complications due to "drying of the respiratory mucosa, such as mucus plugging and endotracheal tube (ETT) occlusion. HMEF is one type of commercial humidification system.
- ◆ Moisture Output : 32mg H<sub>2</sub>O/L
- ◆ Filtration Efficiency : 99.99 %
- ◆ Tidal Volume Range : 150 - 1500 ml
- ◆ Connections : 22M/15F - 22F/15M
- ◆ Gas Monitoring Port : Yes
- ◆ Individually packed in Poly Pack.