


IV Cannula / IV Catheter



- Bevel type needle for easier less traumatic insert.
- Thin Wall PTFE catheter with X- Ray opaque line.
- Flexible and kink resistant catheter is double tapered with rounded tip for smooth introduction and minimum discomfort to the patient.
- Injection port with one way valve and color coded cap for intermittent and safe administration of drugs during infusion.
- Flexible wings provide easy and secure fixation and also prevent skin contact with luer lock connector.
- High transparent flash back chamber for an easy view of returning blood.
- Blood stopper
- Luer lock cap for closing the proximal end of the catheter after removing the needle.

Available Specifications




Code: IT-31

IV Cannula with Luer Lock, Wings and Injection Port
 Sizes: 14G, 16G, 17G, 18G, 20G, 22G & 24G
 Sterile/Blister Pack/Quantities per Carton = 1000 pcs



Code: IT-32

IV Cannula without Wings and Injection Port
 Sizes: 14G, 16G, 17G, 18G, 20G, 22G & 24G
 Sterile/Blister Pack/Quantities per Carton = 1000 pcs



Code: IT-33

IV Cannula with Luer Lock and Wings & without Injection Port
 Sizes: 14G, 16G, 17G, 18G, 20G, 22G & 24G
 Sterile/Blister Pack/Quantities per Carton = 1000 pcs

Size	Color	Catheter Length		Catheter Length (mm)	Water Flow Rate (ml./min.)	Qty. Per Box (pcs.)	Qty. Per Carton (pcs)
		ID mm	OD mm				
14G	Orange	1.7	2.1	45	270	50/100	1000
16G	Grey	1.3	1.7	45	180	50/100	1000
17G	White	1.1	1.5	45	125	50/100	1000
18G	Green	0.9	1.3	45	80	50/100	1000
20G	Pink	0.8	1.1	32	54	50/100	1000
22G	Blue	0.6	0.9	25	33	50/100	1000
24G	Yellow	0.5	0.7	19	20	50/100	1000