



Kehr's T Tube

Code: GE-07

- ◆ Manufactured from ultra soft, flexible, non toxic medical grade PVC for long term use.
- ◆ Sterile, for Aspiration of Mucus and other Secretion from Oropharynx in Newly Born Babies, to ensure Free Respiration.
- ◆ Radio- opaque line provided through out its length for x-ray visualization
- ◆ Sizes : FG: 8, 10, 12, 14, 16,18, 20, 22.

Options Available

- * Various sizes and lengths
- * With or without XRO Line
- * Poly Pack or Blister Pack
- * Inner Bag or Inner Box