

## Laryngeal Mask Airway



Code: AS-27A (PVC)  
Code: AS-27B (Silicon)

- ◆ Cuff of the mask & its tip is anatomically engineered to provide an optimal seal within the hypo-pharynx.
- ◆ Clear transparent tube enable to detect any airway obstruction.
- ◆ Separate inflation line enable the user to move away the line from the teeth of the patient to prevent accidental cuff deflation.
- ◆ Standard 15mm connector facilitates easy connection to the anesthesia circuit.
- ◆ Pilot balloon for easy inflation & deflation of cuff.
- ◆ Made of 100% Silicone or PVC.
- ◆ Available Sizes : 1.0 , 1.5 , 2.0 , 2.5 , 3.0 , 4.0 & 5.0.
- ◆ Individually packed in Paper Poly Pack.