

Nasal Oxygen Catheter



Code: AS-03

- ◆ Designed for administration of oxygen.
- ◆ Green colour catheter made of implantation tested PVC having open distal end with multiple eyes to facilitate dispersion of oxygen and prevent oxygen burn.
- ◆ Colour coded connector for size identification.
- ◆ Sterile, non-toxic & pyrogen free.
- ◆ Available with 8 holes
- ◆ Size (FG) : 6, 8, 10, 12, 14, 16, 18, 20, 22

Oxygen Connecting Tube



Code: AS-13

Options Available

- * Various sizes and lengths
- * With or without XRO Line
- * Poly Pack or Paper-Poly Pack
- * Inner Bag or Inner Box