










- ◆ Made from Non - Toxic Medical Grade PVC Compound
- ◆ Suitable for Gravity Infusions
- ◆ Available in both Vented (for Bottle & Bag) and Non-Vented (for Plastic Bottle) Sharp Spike for penetration.
- ◆ ISO 8536-4 compliant Spike & ISO 594/1 and 594/2 compliant end connection.
- ◆ Available in Latex Free
- ◆ Available in 1/2/3 Needle-Free Injection Site
- ◆ Available in Pre-Slit Injection Site
- ◆ EO Sterile, Non - Toxic and Pyrogen Free

Available Specifications

-  **Code: SP-02-A**
 92" IV Administration Set Non Vented
 Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter - 15micron/Flow Regulator/1 No. Needle-Free Injection 'Y' Site/
 1 No. Latex Free 'Y' Site/Rotating Luer Lock Male Adaptor with Cap/
 2.7mm x 3.9mm PVC Tubing-Length 180+18+18cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 300 pcs
-  **Code: SP-02-B**
 92" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter- 15micron/Flow Regulator/1 No. Needle-Free Injection 'Y' Site/
 1 No. Latex Free 'Y' Site/Rotating Luer Lock Male Adaptor with Cap/
 2.7mm x 3.9mm PVC Tubing-Length 180+18+18cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-C1**
 92" IV Administration Set Non Vented
 Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter-15micron/Flow Regulator/1 No. Needle-Free Injection 'Y' Site/
 Rotating Luer Lock Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing-
 Length 196+18cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-C2**
 100" IV Administration Set Non Vented
 Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter-15micron/Flow Regulator/2 Nos. Needle-Free Injection 'Y' Sites/
 Rotating Luer Lock Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing-
 Length 144+75+15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs



-  **Code: SP-02-C3**
 100" IV Administration Set Non Vented
 Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter - 15micron / Flow Regulator / 2 Nos. Needle - Free
 Injection 'Y' Sites / Back Check Valve / Rotating Luer Lock Male
 Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-Length
 20+124+75+15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-D1**
 83" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter - 15micron/Flow Regulator/Needle-Free
 Injection 'Y' Site/Rotating Luer Lock Male Adaptor with Cap/
 2.7mm x 3.9mm PVC Tubing-Length 180+15cm/
 EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-D2**
 100" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter-15micron/Flow Regulator/2 Nos. Needle-Free
 Injection 'Y' Sites/Back Check Valve/Rotating Luer Lock Male
 Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-
 Length 20+124+75+15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-D3**
 100" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/
 Fluid Filter -15micron/Flow Regulator/3 Nos. Needle-Free
 Injection 'Y' Sites/Rotating Luer Lock Male Adaptor with Cap/
 2.7mm x 3.9mm PVC Tubing-Length 40+95+75+15cm/
 EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-E1**
 83" IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Drip Chamber 60mm/
 Micro Drip Tube/Fluid Filter-15micron/Flow Regulator/
 Needle-Free Injection 'Y' Site/Rotating Luer Lock Male
 Adaptor with Cap/2.7mm x 3.9mm PVC Tubing-
 Length 180+15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-E2**
 100" IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Drip Chamber 60mm/
 Micro Drip Tube/Fluid Filter-15micron/Flow Regulator/
 2 Nos. Needle-Free Injection 'Y' Sites/Back Check Valve/
 Rotating Luer Lock Male Adaptor with Cap/3.0mm x 4.1mm
 PVC Tubing-Length 20+124+75+15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 240 pcs
-  **Code: SP-02-F**
 81" MV Set Vented
 Vented Spike 60 Drops-ml with Handle/Latex Free MV Cylinder
 150ml/Fluid Filter-15micron/Micro Drip Tube/Shut-Off Valve/
 Flow Regulator/Needle-Free Injection 'Y' Site/Slide Clamp
 Rotating Luer Lock Male Adaptor with Cap/2.7mm x 3.9mm
 PVC Tubing-Length 25+135+15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 80 pcs