

PACLITAXEL ADMINISTRATION SETS

- ◆ Self and fast priming of adult I.V. Set to minimize nursing time.
- ◆ Automatic air-venting in any position loss of flow due to air locks.
- ◆ Able to withstand a minimum of 45psi (3.1 bar) pressure.
- ◆ ISO 8536-4 compliant Spike & ISO 594/1 and 594/2 compliant end connection.
- ◆ Bacterial retention throughout 96 hrs. Flow.
- ◆ Options available with PVC/DEHP Free
- ◆ Options available with Filter media 1.2 / 0.2 micron
- ◆ EO Sterile, Non - Toxic and Pyrogen Free



Infusion Therapy

Available Specifications



Code: SP-04

76" Infusion for Paclitaxel Administration
Vented Spike 20 drops-ml with Guard/PVC-DEHP Free Drip Chamber 60mm/
Fluid Filter - 15micron/Flow Regulator/0.2 Micron Filter/Rotating Luer
Lock Male Adaptor with Cap/3.0mm x4.1mm PVC-DEHP Free Tubing-
Length 140+15cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 300 pcs

Also available with Latex Free 'Y' Site

PACLITAXEL Administration