


IV Set - Standard Models (Paediatric)




Code: IT-26

63" Paediatric IV Administration Set Non-Vented
 Non-Vented Spike 60 drops-ml with Guard/Micro Drip Tube/Drip Chamber 50mm with /
 Fluid Filter/Flow Regulator/Tube Latex/Luer Tip Male Adaptor with Cap/Needle 23g x 1"
 with Protector/2.8mm x 3.8mm PVC Tubing- Length 150cm/
 EO Sterile/Poly Pack or Paper Poly Pack
 Quantities per Carton = 600 pcs in Poly Pack or 400 pcs in Paper Poly Pack




Code: IT-38

65" Paediatric IV Administration Set Non-Vented
 Non-Vented Spike 60 drops-ml with Guard/Micro Drip Tube/Drip Chamber 60mm/Fluid Filter/
 Latex Free 'Y' Site/Flow Regulator/Luer Lock Male Adaptor with Cap/Needle 23g x 1"
 with Protector/2.7mm x 3.9mm PVC Tubing- Length 135+15cm/EO Sterile/ Paper Poly Pack
 Quantities per Carton = 350 pcs



Code: IT-39

63" Paediatric IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Micro Drip Tube/Drip Chamber 50mm/Fluid Filter/
 Flash Bulb/Flow Regulator/Luer Tip Male Adaptor with Cap/Needle 23g x 1" with Protector
 2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack or Paper Poly Pack
 Quantities per Carton = 560 pcs in Poly Pack or 400 pcs in Paper Poly Pack



Code: IT-40

65" Paediatric IV Administration Set Vented
 Vented Spike with 60 drops-ml Guard/Drip Chamber 60mm/Micro Drip Tube/Fluid Filter/
 Latex Free 'Y' Site/Flow Regulator/Luer Lock Male Adaptor with Cap/Needle 23g x 1"
 with Protector/2.7mm x 3.9mm PVC Tubing- Length 135+15cm/EO Sterile/ Paper Poly Pack
 Quantities per Carton = 350 pcs



Code: IT-42

65" Paediatric IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Drip Chamber 50mm with Micro Drip Tube/Fluid
 Filter/Latex Free 'Y' Injection Site/Flow Regulator/Luer Lock Male Fitting/Needle 23G with
 Protector/2.8mm x 3.8mm PVC Tubing-Length 135+15cm/EO Sterile/Poly Pouch
 Quantities Per Carton = 560 pcs



Code: IT-43

65" Paediatric IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Drip Chamber 50mm with Micro Drip Tube/Fluid
 Filter/Bulb Latex/Flow Regulator/Luer Lock Male Fitting/Needle 23G with Protector/
 2.8mm x 3.8mm PVC Tubing-Length 150cm/EO Sterile/Poly Pouch
 Quantities Per Carton = 560 pcs