

Peritoneal Dialysis Catheter Kit



Code: IT-29A

- ◆ Stainless steel trocar is provided with triple angle sharpness to facilitate smooth penetration.
- ◆ Spare scalpel is supplied to provide the incision.
- ◆ Multi hole, open end catheter is manufactured from medical grade polyamide.
- ◆ The trocar and the tip of the catheter is perfectly matched to facilitate trauma free entry.
- ◆ Latex free flash ball is provided on the junction unit for additional medication.
- ◆ Unit consists:
 - *Stainless steel trocar
 - *Multi hole smooth tip catheter
 - *Junction unit with control clamp
 - *Scalpel
- ◆ Available Size : Adult & Child