




## Peritoneal Dialysis Transfusion Set

Infusion / Transfusion Therapy

- Manufactured from Non-Toxic & Non-Irritant PVC.
- Extra soft kink-resistant frosted tubing.
- Specially designed 'Y' shape transfusion set for administering dialysis solutions.
- Provided with two upper control clamps to facilitate the alternative change of solution bottles.
- Two lower control clamps help easy input and drainage of solution.
- Proximal end is provided with suitable adapter for easy connection to catheter unit.
- Sterile Pack - ETO (Ethylene Oxide) Sterilization validate according to Harmonized Standards



Code: IT-29

Peritoneal Dialysis Set