

Infusion Set for Photosensitive Drugs

- Made Lightproof/Photosensitive polymer materials
- Available in Yellow and Brown colour
- Available with or without Precision Flow Controller [Dial type] 5 ml to 250 ml / hrs.
- Latex Free
- Soft Clear Flexible Tubing
- Sterile, Non-Toxic and Pyrogen-Free.
- Sterilized by E.O. (Ethylene Oxide)

Available Specifications



Code: SP-08A

88" IV Administration Set Vented

Vented Spike 20 Drops-ml with Guard/Brown Colour Drip Chamber 60mm/Fluid Filter 15 Micron/
Flow Regulator/Precision Flow Controller 5-250ml-hr./Latex Free 'Y' Site/Rotating Luer Lock

Male Adaptor with Cap/3.0 x 4.1mm PVC Tubing Brown Colour -

Length 30+140+30cm/Sterile/Paper Poly Pack/

Quantities per Carton = 240 pcs



Code: SP-08B

88" IV Administration Set Vented

Vented Spike 20 Drops-ml with Guard/Yellow Colour Drip Chamber 60mm/Fluid Filter 15 Micron/
Flow Regulator/Latex Free Flash Bulb/Rotating Luer Lock Male Adaptor with Cap/

3.0 x 4.1mm PVC Tubing Yellow Colour - Length 150cm/Sterile/Paper Poly Pack/

Quantities per Carton = 300 pcs

