




## Infusion Sets - PVC/DEHP Free

- Pressure resistant Infusion Set
- PVC / DEHP Free Tubing and Chamber
- Sharp spike for better penetration.
- ISO 8536-4 compliant spike ISO 594/1 & ISO 594/2 compliant end-connection
- 15 Micron Fluid Filter
- Non-toxic, Pyrogen free
- Gamma / ETO Sterile

### Available Specifications



Code: SP-06-A

60" Pressure Resistant IV Administration Set Vented  
 Vented Spike 20 Drops-ml with Guard/Drip Chamber 60mm/  
 Fluid Filter -15micron/Flow Regulator/Male Luer Lock Adaptor  
 with Cap/Back Check Valve/3.0mm x 4.1mm Tubing-  
 Length 150cm/EO Sterile/Paper Poly Pack  
 Quantities per Carton = 300 pcs



Code: SP-06-B

60" IV Administration set Vented  
 Vented Spike 20 Drops-ml with Guard/Drip Chamber 60mm/  
 Fluid Filter - 15micron/Flow Regulator/Male Luer Lock  
 Adaptor with Cap/3.0mm x 4.1mm Tubing-Length 150cm/  
 EO Sterile/Paper Poly Pack  
 Quantities per Carton = 300 pcs



Code: SP-06-C

70" IV Administration IV Set Vented  
 Vented Spike 20 Drops-ml with Guard/Drip Chamber 60mm/  
 Fluid Filter - 15micron/Latex Free 'Y' Site/Flow Regulator/  
 Male Luer Lock Adaptor with Cap/3.0mm x 4.1mm Tubing-  
 Length 150+15cm/EO Sterile/Paper Poly Pack  
 Quantities per Carton = 300 pcs



Code: SP-06-D

60" IV Administration Set Vented  
 Vented Spike 20 Drops-ml with Guard/Drip Chamber 60mm/  
 Fluid Filter - 15micron/Latex Free Flash Bulb/Flow Regulator/  
 Male Luer Lock Adaptor with Cap/3.0mm x 4.1mm Tubing-  
 Length 150cm/EO Sterile/Paper Poly Pack  
 Quantities per Carton = 300 pcs

Infusion Therapy

PVC/DEHP FREE