

## Ryle's Tube



Code: GE-02

Specially designed for patient undergoing long term gastro intestinal feeding, aspiration and intestinal secretion.

Siliconized and smooth outer surface to reduce chances of tissue irritation. Thermosensitive material softening with body temperature conforming to the shape of naso - gastric tract.

Corrosion resistant stainless steel balls sealed at the distal end of tube for X-ray identification and for assisting the movement of tube during insertion.

Provided with full length radio - opaque line for checking the tube placement.

Provision of four lateral eyes for efficient aspiration.

Marked with 1, 2, & 3 black circles at 500mm, 600mm & 700mm for accurate placement

Provided with soft tapered connector at the proximal end to fit easily with standard accessories.

Sterile by Ethylene Oxide

Sizes(FG) : 8, 10, 12, 14, 16, 18, 20, 22, 24, 26

## Deodenal Ryle's Tube



Code: GE-02-A

Nasogastric introduction for aspiration of intestinal fluid or for administration of enteral nutrition.

Large lateral eyes permit effective drainage

Sizes of eyes compatible with lumen preventing blockage

Kink-resistant

Transparent material facilitates observation of aspirate

Frosted surface eliminates the need for lubricant

X-ray line and tip (ryles) ensure position control

Colour coded funnels with female connector and closure fit standard adaptors

Length 105cm, Size: FG: 8 - 20

Tip X-ray detectable during passage to duodenum

PVC, Medical grade

Single packed, Sterile by Ethylene Oxide

### Options Available

- \* Various sizes and lengths
- \* With or without XRO Line
- \* Poly Pack or Blister Pack
- \* Inner Bag or Inner Box