


Infusion Sets For Secondary Medication


- Ideally suited for Secondary Medication infusions
- Available in both vented (for bottles and bags) and non-vented (for Plastic Bottle bottles) Sharp spike for better penetration
- ISO 8536-4 compliant spike ISO 594/1 & ISO 594/2 compliant end-connection
- Soft Clear Flexible Tubing
- Sterile, Non-Toxic and Pyrogen-Free.
- Sterilized by E.O. (Ethylene Oxide) Gas



Infusion Therapy


For Secondary Medication

Available Specifications




Code: SP-07-A

32" IV Administration Set Vented
 Vented Spike 20 Drops-ml with Guard/Drip Chamber 60mm/Fluid Filter - 15micron/
 Flow Regulator/Male Luer Tip Adaptor with Cap/3.0mm x 4.1mm Tubing-
 Length 70cm/Hanger /EO Sterile/Paper Poly Pack
 Quantities per Carton = 350 pcs



Code: SP-07-B

32" IV Administration Set Non-Vented
 Non-Vented Spike 20 Drops-ml with Guard/Drip Chamber 60mm/Fluid Filter- 15micron/
 Flow Regulator/Male Luer Tip Adaptor/Needle 18g with Protector/
 3.0mm x 4.1mm Tubing-Length 70cm/Hanger /EO Sterile/Paper Poly Pack
 Quantities per Carton = 350 pcs



Code: SP-07-C

32" IV Administration Set Non-Vented
 Non-Vented Spike 10 Drops-ml with Guard/Drip Chamber 60mm/Fluid Filter -
 15micron/Flow Regulator/Male Luer Tip Adaptor with Cap/3.0mm x 4.1mm
 Tubing-Length 70cm/Hanger /EO Sterile/Paper Poly Pack
 Quantities per Carton = 350 pcs