


BT Set - Standard Models



- Sharp spike vented (for bottles & bags) and non-vented (for plastic bottles) for better penetration.
- ISO 8536-4 compliant spike, ISO 594/1 & ISO 594/2 compliant end-connection
- Superior quality of latex tubing/flash bulb for self-sealing.
- 200 micron filter and standard size surface area (10 sq. Cm)
- Specially designed roller clamp for better flow control.
- Top quality needle.
- Available with or without air-vent
- Non-Toxic, Pyrogen Free & Sterile by EO

Transfusion Therapy - Standard Models

Available Specifications



Code: TT-03

64" Blood Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard-Bevel Type/Drip Chamber 100mm/Blood Filter 60mm/Bulb Latex/Flow Regulator/Luer Tip Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack or Paper Poly
 Quantities per Carton = 400 pcs in Poly Pack & 325 pcs in Paper Poly



Code: TT-04

64" Blood Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 100mm/Blood Filter 60mm/Bulb Latex/Flow Regulator/Luer Tip Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 400 pcs



Code: TT-05

64" Blood Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 100mm/Blood Filter 60mm/Latex Free 'Y' Site/Flow Regulator/Luer Lock Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 300 pcs



Code: TT-06

64" Blood Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard-Bevel Type/Drip Chamber 100mm/Blood Filter 60mm/Bulb Latex/Flow Regulator/Luer Lock Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 325 pcs



Code: TT-07

64" Blood Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard-Bevel Type/Drip Chamber 100mm/Blood Filter 60mm/Latex Free 'Y' Site/Flow Regulator/Rotating Luer Lock Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 300 pcs

BT Set - Standard Models




Code: TT-08

64" Blood Administration Set Non-Vented
 Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 100mm/Blood Filter 60mm/
 Latex Free 'Y' Site/Flow Regulator/Rotating Luer Lock Adaptor with Cap/2.7mm x 3.9mm PVC
 Tubing- Length 150cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 300 pcs



Code: TT-09

77" Blood Administration Set Non-Vented
 Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 100mm + 70mm (Double Chamber)/
 Blood Filter 60mm/Latex Free 'Y' Site/Flow Regulator/Rotating Luer Lock Adaptor with Cap/
 2.7mm x 3.9mm PVC Tubing- Length 165 + 15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 175 pcs




Code: TT-10

64" Blood Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Tube Latex /
 Flow Regulator/ Luer Tip Male Adaptor/Needle 18g x 1.5" with Protector/
 2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 500 pcs



Code: TT-11

64" Blood Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Tube Latex /
 Flow Regulator/ Luer Tip Male Adaptor/Needle 18g x 1.5" with Protector/
 2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 500 pcs




Code: TT-13

64" Blood Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Bulb Latex /
 Flow Regulator/ Luer Tip Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
 Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 500 pcs



Code: TT-14

64" Blood Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Bulb Latex /
 Flow Regulator/ Luer Tip Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
 Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 450 pcs



Code: TT-15

101" 'Y' Division Blood Administration Set
 2 Nos. Non-Vented Spike 15 drops-ml with Guard-Bevel type/2 Nos. Blood Chamber 100mm/
 1 no. Blood Filter 60mm/Latex Free 'Y' Site/PVC 2Way Connector/3 Nos. Flow Regulators/
 Rotating Luer Lock Male Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-
 Length 2 nos. 50+135+15cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 150 pcs