

IV Set - Standard Models



- Sharp spike vented (for bottles & bags) and non-vented (for plastic bottles) for better penetration.
- ISO 8536-4 compliant spike, ISO 594/1 & ISO 594/2 compliant end-connection
- Superior quality of latex tube/flash bulb for self-sealing.
- Specially designed roller clamp for better flow control.
- Top quality needle.
- Available with or without air-vent
- Option available with Latex Free
- Non-Toxic, Pyrogen Free & Sterile by EO

Infusion Therapy - Standard Models

Available Specifications



Code: IT-04

63" IV Administration Set Vented

Vented Spike-20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/Tube Latex/
Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.25" with Protector/
2.7mm x 3.6mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 600 pcs



Code: IT-05

63" IV Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/Bulb Latex/
Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.25" with Protector/
2.7mm x 3.6mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 560 pcs



Code: IT-05L

63" IV Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/Bulb Latex/
Big Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.25" with Protector/
2.7mm x 3.6mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 560 pcs



Code: IT-06

62" IV Administration Set Non-Vented with Butterfly Needle

Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/
Tube Latex/Flow Regulator/Luer Tip Male Adaptor/12" x 21g Scalp Vein Set/
2.8mm x 3.8mm PVC Tubing- Length 120cm/EO Sterile/Poly Pack
Quantities per Carton = 600 pcs




Code: IT-07

63" IV Administration Set Vented


Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/Latex Free
'Y' Site/Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.5" with Protector/
2.8mm x 3.8mm PVC Tubing- Length 135+15cm/EO Sterile/Poly Pack
Quantities per Carton = 560 pcs

IV Set - Standard Models



Code: IT-08

63" IV Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/
 Flash Bulb/Flow Regulator/Luer Lock Male Adaptor with Cap/2.8mm x 3.8mm
 PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 560 pcs




Code: IT-09

63" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/Flash Bulb/
 Big Flow Regulator/Luer Lock Male Adaptor/Needle 21g x 1.5" with Protector/
 2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Paper Poly Pack
 Quantities per Carton = 350 pcs




Code: IT-09A

63" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/Latex Free 'Y' Site
 Big Flow Regulator/Luer Lock Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
 Length 135+15cm/EO Sterile/ Poly Pack
 Quantities per Carton = 560 pcs




Code: IT-10

63" Bottle Needle IV Administration Set Non-Vented
 16g x 1.25" Bottle Needle Single Slot with Cap 20 drops-ml/Drip Chamber 50mm/
 Fluid Filter/Tube Latex/Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.5"
 with Protector/3.0mm x 4.1mm PVC Tubing- Length 15cm/2.8mm x 3.8mm
 PVC Tubing- Length 135cm/Extra Air-Way Needle 16g x 1.25" Single Slot
 Needle with Cotton and 5cm long PVC Tubing/EO Sterile/Poly Pack
 Quantities per Carton = 600 pcs



Code: IT-11

63" IV Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/
 Latex Free 'Y' Site/Luer Tip Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
 Length 135+15cm/Flow Regulator/EO Sterile/Poly Pack or Paper Poly Pack
 Quantities per Carton = 560 pcs in Poly or 400 pcs in Paper Poly



Code: IT-12

63" IV Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/
 Tube Latex/Luer Tip Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
 Length 150cm/Flow Regulator/EO Sterile/Poly Pack or Paper Poly Pack
 Quantities per Carton = 560 pcs in Poly or 400 pcs in Paper Poly



Code: IT-13

63" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/
 Latex Free 'Y' Site/Luer Lock Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
 Length 135+15cm/Flow Regulator/EO Sterile/Paper Poly Pack
 Quantities per Carton = 400 pcs



Code: IT-14

65" IV Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/
 Flash Bulb/Luer Tip Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing-
 Length 150cm/Big Flow Regulator/EO Sterile/Paper Poly Pack
 Quantities per Carton = 400 pcs

IV Set - Standard Models



Code: IT-15

72" IV Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/Latex Free 'Y' Site/Big Flow Regulator/Rotating Luer Lock Male Adaptor with Cap/ 2.7mm x 3.9mm PVC Tubing- Length 150+15cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 350 pcs



Code: IT-15-L

72" IV Administration Set Vented


Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter 15 micron/ Latex Free 'Y' Injection Site/Big Flow Regulator/Luer Lock Male Fitting/Needle 21G x 1.5" with Protector/2.7mm x 3.9mm PVC Tubing-Length 150+15cm/EO Sterile/Poly Pouch
Quantities Per Carton = 500 pcs



Code: IT-16

63" IV Administration Set Vented


Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/Latex Free 'Y' Site/Luer Lock Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing- Length 135+15cm/Flow Regulator/EO Sterile/Poly Pack or Paper Poly Pack
Quantities per Carton = 560 pcs in Poly Pack or 400 pcs in Paper Poly Pack



Code: IT-17

78" I.V. Administration Set Vented


Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/Latex Free 'Y' Site - 2 nos./Big Flow Regulator/Rotating Luer Lock Male Adaptor with Cap/ 2.7mm x 3.9mm PVC Tubing- Length 50+100+30cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 300 pcs



Code: IT-18

66" IV Administration Set Vented with 3 Way Stopcock

Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/ Latex Free 'Y' Site/Big Flow Regulator/2 nos. Rotating Luer Lock Male Adaptor with Cap/3 Way Stop Cock/2.7mm x 3.9mm PVC Tubing- Length 150+15+15cm/ EO Sterile/Paper Poly Pack
Quantities per Carton = 200 pcs



Code: IT-19

78" IV Administration Set Vented with Back Check Valve

Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/Back Check Valve/2 Nos. Latex Free 'Y' Site/Flow Regulator/Rotating Luer Lock Male Adaptor with Cap/ 2.7mm x 3.9mm PVC Tubing- Length 25+25+100+30cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 200 pcs



Code: IT-21

63" IV Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/ Flash Bulb/Flow Regulator/Luer Tip Male Adaptor /Needle 21g x 1.5" with Protector/ 2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 500 pcs

IV Set - Standard Models



Code: IT-21L

63" IV Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 50mm/Fluid Filter/
Flash Bulb/Big Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.5" with
Protector/ 2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack

Quantities per Carton = 500 pcs




Code: IT-22

65" IV Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/
Latex Free 'Y' Site/Big Flow Regulator/Rotating Luer Lock Male Adaptor with Cap/
2.7mm x 3.9mm PVC Tubing- Length 135+15cm/EO Sterile/Paper Poly Pack

Quantities per Carton = 350 pcs




Code: IT-23

65" IV Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/
Big Flow Regulator/Latex Free 'Y' Site/Rotating Luer Lock Male Adaptor with Cap/
2.7mm x 3.9mm PVC Tubing- Length 135+15cm/EO Sterile/Paper Poly Pack

Quantities per Carton = 300 pcs



Code: IT-24

65" IV Administration Set Vented


Vented Spike 20 drops-ml with Guard/Drip Chamber 60mm/Fluid Filter/Latex Free
'Y' Site with Wing/Big Flow Regulator/Rotating Luer Lock Male Adaptor with Cap/
2.7mm x 3.9mm PVC Tubing- Length 135+15cm/EO Sterile/Paper Poly Pack

Quantities per Carton = 300 pcs

Options

- **Needle-Free**
- **PVC/DEHP Free**
- **With Check Valve**
- **For Infusion Pump**
- **For Filter Administration**
- **For Photosensitive Drugs**
- **For Secondary Medication**
- **For Oncology/Taxol Solutions**

IV Set - Standard Models (Paediatric)



Code: IT-26

63" Paediatric IV Administration Set Non-Vented
 Non-Vented Spike 60 drops-ml with Guard/Micro Drip Tube/Drip Chamber 50mm /
 Fluid Filter/Flow Regulator/Tube Latex /Luer Tip Male Adaptor/Needle 23g x 1"
 with Protector/2.8mm x 3.8mm PVC Tubing- Length 150cm/
 EO Sterile/Poly Pack or Paper Poly Pack
 Quantities per Carton = 600 pcs in Poly Pack or 400 pcs in Paper Poly Pack



Code: IT-38

65" Paediatric IV Administration Set Non-Vented
 Non-Vented Spike 60 drops-ml with Guard/Micro Drip Tube/Drip Chamber 60mm/Fluid Filter/
 Latex Free 'Y' Site/Big Flow Regulator/Luer Lock Male Adaptor/Needle 23g x 1"
 with Protector/2.7mm x 3.9mm PVC Tubing- Length 135+15cm/EO Sterile/ Paper Poly Pack
 Quantities per Carton = 350 pcs




Code: IT-39

63" Paediatric IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Micro Drip Tube/Drip Chamber 50mm/Fluid Filter/
 Flash Bulb/Flow Regulator/Luer Tip Male Adaptor/Needle 23g x 1" with Protector
 2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack or Paper Poly Pack
 Quantities per Carton = 560 pcs in Poly Pack or 400 pcs in Paper Poly Pack



Code: IT-40

65" Paediatric IV Administration Set Vented
 Vented Spike 60 drops-ml Guard/Drip Chamber 60mm/Micro Drip Tube/Fluid Filter/
 Latex Free 'Y' Site/Big Flow Regulator/Luer Lock Male Adaptor/Needle 23g x 1"
 with Protector/2.7mm x 3.9mm PVC Tubing- Length 135+15cm/EO Sterile/ Paper Poly Pack
 Quantities per Carton = 350 pcs



Code: IT-42

65" Paediatric IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Drip Chamber 50mm with Micro Drip Tube/Fluid
 Filter/Latex Free 'Y' Injection Site/Flow Regulator/Luer Lock Male Fitting/Needle 23Gx1" with
 Protector/2.8mm x 3.8mm PVC Tubing-Length 135+15cm/EO Sterile/Poly Pouch
 Quantities Per Carton = 560 pcs



Code: IT-43

65" Paediatric IV Administration Set Vented
 Vented Spike 60 drops-ml with Guard/Drip Chamber 50mm with Micro Drip Tube/Fluid
 Filter/Bulb Latex/Flow Regulator/Luer Lock Male Fitting/Needle 23Gx1" with Protector/
 2.8mm x 3.8mm PVC Tubing-Length 150cm/EO Sterile/Poly Pouch
 Quantities Per Carton = 560 pcs