







## Scalp Vein Set

- Flexible fixing wings (butterfly) needle for pediatric and geriatric use.
- Soft wing for less pain to patient
- Ideal for single push injections
- Gauge sizes available include: 20, 21, 22, 23, 24 & 25
- Siliconized, ultra sharp, stainless steel, bevel tip needle
- Effective needle length : 20 mm
- Color coded
- 30 cm, kink resistant tubing
- Luer tip female connector with closure

### Available Specifications

20G		Yellow
21G		Green
22G		Black
23G		Blue
24G		Purple
25G		Orange



Code: IT-28

