



Code: UR-05

## 'Y' Junction Set

'Y' Junction Set / T.U.R. Set (Trans Urethral Resection Set) designed for continuous irrigation during trans urethral resection.

Y shaped connector facilitates alternative changes of solutions.

Slide control clamp aids quick changes of bottle.

Proximal end fitted with flexible latex tube allows easy connection to endoscope.

E.O. Sterile, Non-toxic and Pyrogen-free.

### Options Available

- \* Various sizes and lengths
- \* With or without XRO Line
- \* Poly Pack or Paper-Poly Pack
- \* Inner Bag or Inner Box