

Urine Pot (2 in 1)



Code: UR-50

Manufacturing is carved out of the best quality raw material. These can hygienically be used in hospital wards and various other patient areas for urine collection from a patient necessitating the use of urine pot.

- 2 in 1 Urine Pot for Male & Female also.
- 1000 ml Capacity with graduation mark,
- Easy fit cover for Hygiene and Odourless,
- Easy carry handle.
- Can be stored in standing or lying position.
- Disposable. Available in loose pack/ shipper pack.
- Individual Packed. / Loose pack. & Shipper pack also.
- Packing Details: 595 X 375 X 325 mm.-→ 25 Nos.