

# Urine Meter

## Three Chambers System

Measuring glass: 400ml  
 Urine collection bag for convenient drainage of urine.  
 Inlet tube 120cm-170cm  
 Sampling Port  
 Bed-Sheet Clamp  
 Drainage control knob effectively drains urine.  
 Hanger & string for convenient and firm fixing of the urine meter to the upper side frame of the bed.  
 Anti-reflux drip chamber.  
 Exact and fine scaling for timed measurements  
 Urine sampling points at connector and measuring glass  
 Germs- and liquid tight, hydrophobic coated filter for pressure equalisation  
 Long, highly flexible and bend-protected tube  
 Single hand, easily operated outlet valve  
 Sterile by EO

**Latex Free**



Code: UR-12

## Single Chamber System



Code: UR-13

Measuring glass: 500ml  
 Urine collection bag for convenient drainage of urine.  
 Inlet tube 120cm-170cm  
 Sampling Port  
 Bed-Sheet Clamp  
 Drainage control knob effectively drains urine.  
 Hanger & string for convenient and firm fixing of the urine meter to the upper side frame of the bed.  
 Anti-reflux drip chamber.  
 Exact and fine scaling for timed measurements  
 Urine sampling points at connector and measuring glass  
 Germs- and liquid tight, hydrophobic coated filter for pressure equalisation  
 Long, highly flexible and bend-protected tube  
 Single hand, easily operated outlet valve  
 Sterile by EO

**Latex Free**

## Single Chamber System

Measuring glass: 200 / 250ml  
 Urine collection bag for convenient drainage of urine.  
 Inlet tube 120cm-170cm  
 Drainage control knob effectively drains urine.  
 Hanger & string for convenient and firm fixing of the urine meter to the upper side frame of the bed.  
 Exact and fine scaling for timed measurements  
 Germs- and liquid tight.  
 Highly flexible and bend-protected tube  
 Single hand, easily operated outlet valve  
 Sterile EO

**Latex Free**



Code: UR-09