



Nelaton Catheter

Code: UR-01

Designed for short term bladder catheterisation.

Made from medical grade thermosensitive PVC material. Non-irritant to delicate mucous membranes of urethra.

Perfectly finished closed distal end for smooth & painless insertion.

Provided with two side eyes for efficient drainage. Provided with funnel shaped color-coded connector at the proximal end.

With X-ray opaque line made out of tested and siliconised, non-toxic PVC material with smooth finish to prevent obstruction of lumen and trauma to urethra by abrasion.

Sterile, non-toxic and pyrogen free.

Size (FG): 8, 10, 12, 14, 16, 18, 20, 22, 24

Options Available

- * Various sizes and lengths
- * With or without XRO Line
- * Poly Pack or Paper-Poly Pack
- * Inner Bag or Inner Box