



Code: AS-14

Adjustable Venturi Mask

Made of clear, Non-Toxic PVC
Colour coded kit for ease of use and delivery
of variable concentrations of oxygen.
Two Color coded diluters for seven concentrations:
Green diluter for 24%, 26%, 28%, 30%
White diluter for 35%, 40%, 50%
Six Color coded diluters for six concentrations:
24%, 28%, 31%, 35%, 40%, 50%
Adjustable nose clip
Adaptor for high humidity
With 2m anti-crush tube.

Options Available

- * Adult & Child
- * Poly Pack or Paper-Poly Pack
- * Inner Bag or Inner Box



Code: AS-14