

## Endotracheal Tube Plain

Code: AS-08

Pilot balloon with integral one way valve

Size of the tube clearly marked on the tube

Cuff is optimally designed for routine intubation


One way valve for efficient and easy inflation and deflation of cuff

Murphy eye, high volume, low pressure cuff, oral/nasal

ETO Sterile, Latex free

Available sizes: 2.0mm to 10.0mm

Paper Poly Packing



## Endotracheal Tube Cuffed

Code: AS-07

Non-return valve with indicator facilitate the cuff inflation and deflation

Thin, smooth cuff inflates symmetrically eliminating any risk of herniation

Terminal end is suitable for 15mm connector

Individual Packing in Blister Pack

ETO Sterile, Latex free

Available sizes: 3.0mm to 10mm

Paper Poly Packing