



Code: TT-30AV

AV FISTULA

- ◆ AV Fistula Needle allows Arterial and Venous access for patients undergoing haemodialysis
- ◆ Flexible rotating wings designed for tight gripping during intuition and easy location, and coded for size identification
- ◆ Needle with bevel edge with automated siliconization
- ◆ Minimizes pain during insertion, maximizes smooth, laminar flow
- ◆ UPS Class VI compatible kink resistant tubing
- ◆ Luer has universal 6% taper allow smooth connection with standard bloodline sets
- ◆ Available in single pack (A/V) or twin pack (A&V)
- ◆ Available with or without Back-Eye Needle
- ◆ Available with Red, Blue or White Clamp
- ◆ Available with 15cm or 30cm Tubing
- ◆ Sizes available: 15G, 16G & 17G

