


Breathing Circuit with Single & Double Water Trap



Code: AS-25A (Single Water Trap)
Code: AS-25B (Double Water Trap)

- ◆ Re-sealable water-trap provides a perfect sealing while removal of excess collection of water from the water-traps giving a continues flow.
- ◆ Made of high quality medical grade plastic with 22M-22F Straight Connector.
- ◆ Transparent for easy monitoring of water build-up or blockages.
- ◆ Corrugated, Flexible and Lightweight tubing 1.6 mtr. with single & double re-sealable Water Trap, 'Y piece' with monitoring port.
- ◆ Available Sizes : Adult & Pediatric .
- ◆ Individually packed in Poly Pack .