


Double Swivel Mount (Catheter T Mount)



Code: AS-29

- Double swivel mount connector with suction port and cap provides better flexibility, easy accessibility and provide comfort to patient during suction. It is also known as Catheter T Mount.
- 15 mm I.D./22 mm O.D. connector fits all tracheal tubes.
- 22 mm standard female connector is compatible with all breathing & ventilator circuits.
- Specially designed port at the elbow facilitates suction & gas sampling without disconnecting the breathing circuit.
- Collapsible corrugated tubing is inert to all anesthetic gases & reagents.
- Individually packed in Paper Poly Pack.