


IV Set - Economy Models



- Sharp spike vented (for bottles & bags) and non-vented (for plastic bottles) for better penetration.
- ISO 8536-4 compliant spike, ISO 594/1 & ISO 594/2 compliant end-connection
- Superior quality of latex tube/flash bulb for self-sealing.
- Specially designed roller clamp for better flow control.
- Top quality needle.
- Available with or without air-vent
- Non-Toxic, Pyrogen Free & Sterile by EO

Infusion Therapy - Economy Models

Available Specifications



Code: IT-01

63" IV Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 45mm/Fluid Filter/Tube Latex/
 Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.25" with Protector/
 2.6mm x 3.2mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 640 pcs



Code: IT-01A

63" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 45mm/Fluid Filter/Tube Latex/
 Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.25" with Protector/
 2.6mm x 3.2mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 640 pcs



Code: IT-05C

63" IV Administration Set Vented
 Vented Spike 20 drops-ml with Guard/Drip Chamber 45mm/Fluid Filter/Bulb Latex/
 Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.25" with Protector/
 2.6mm x 3.2mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 640 pcs



Code: IT-03

63" IV Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 45mm/Fluid Filter/
 Tube Latex/Flow Regulator/Luer Tip Male Adaptor/Needle 21g x 1.25" with Protector/
 2.7mm x 3.6mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
 Quantities per Carton = 640 pcs



Code: IT-41

57" IV Administration Set Non-Vented
 Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 45mm/Fluid Filter/
 Flow Regulator/Luer Tip Male Adaptor with Cap/2.7mm x 3.6mm
 PVC Tubing- Length 135cm/EO Sterile/Poly Pack
 Quantities per Carton = 640 pcs