



Code: GE-03

Levin's Tube

- ◆ Designed for nasal & gastric introduction for nutrition & aspiration of intestinal secretion.
- ◆ Siliconised & smooth outer surface to reduce chances of tissue irritation. Thermosensitive material softening with body temperature conforming to the shape of naso-gastric tract.
- ◆ Provided with full-length radio-opaque line for checking the tube placement.
- ◆ Provision of three eye-holes for efficient aspiration.
- ◆ Soft tapered connector at the proximal end to fit easily with standard accessories.
- ◆ Sterile, non-toxic and pyrogen free.
- ◆ Size (FG): 8, 10, 12, 14, 16, 18, 20, 22, 24

Options Available

- * Various sizes and lengths
- * With or without XRO Line
- * Poly Pack or Blister Pack
- * Inner Bag or Inner Box