

Nebulizer Mask Kit with 'T' Mouth



Code : AS-31

- ◆ For delivery of medication through 'T' Mouthpiece
- ◆ T-piece may be sited between the Y-piece of the breathing system and Tracheal tube / Catheter mount
- ◆ Threaded, smooth Locking system to prevent Leakages
- ◆ Consistent performance in both Vertical and Sleeping position
- ◆ 'CC' marking for medicine measurement
- ◆ Star Lu-men Crush Resistant Oxygen Tubing - 2 mtr
- ◆ Optimal particle size delivered at efficient Nebulization rate