

BT Set - Standard Models



- Sharp spike vented (for bottles & bags) and non-vented (for plastic bottles) for better penetration.
- ISO 8536-4 compliant spike, ISO 594/1 & ISO 594/2 compliant end-connection
- Superior quality of latex tubing/flash bulb for self-sealing.
- 200 micron filter and standard size surface area (10 sq. Cm)
- Specially designed roller clamp for better flow control.
- Top quality needle.
- Available with or without air-vent
- Non-Toxic, Pyrogen Free & Sterile by EO

Transfusion Therapy - Standard Models

Available Specifications



Code: TT-03

64" Blood Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard-Bevel Type/Drip Chamber 100mm/Blood Filter 60mm/Bulb Latex/Flow Regulator/Luer Tip Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack or Paper Poly
Quantities per Carton = 400 pcs in Poly Pack & 325 pcs in Paper Poly

Option: (TT-03L) Above Specification with Needle & Poly Pack



Code: TT-04

64" Blood Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 100mm/Blood Filter 60mm/Bulb Latex/Flow Regulator/Luer Tip Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 400 pcs



Code: TT-05

64" Blood Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 100mm/Blood Filter 60mm/Latex Free 'Y' Site/Flow Regulator/Luer Lock Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 300 pcs



Code: TT-06

64" Blood Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard-Bevel Type/Drip Chamber 100mm/Blood Filter 60mm/Bulb Latex/Flow Regulator/Luer Lock Male Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 325 pcs



Code: TT-07

64" Blood Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard-Bevel Type/Drip Chamber 100mm/Blood Filter 60mm/Latex Free 'Y' Site/Flow Regulator/Rotating Luer Lock Adaptor with Cap/2.7mm x 3.9mm PVC Tubing- Length 150cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 300 pcs

BT Set - Standard Models



Code: TT-08

64" Blood Administration Set Non-Vented

Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 100mm/Blood Filter 60mm/
Latex Free 'Y' Site/Flow Regulator/Rotating Luer Lock Adaptor with Cap/2.7mm x 3.9mm PVC
Tubing- Length 150cm/EO Sterile/Paper Poly Pack / Quantities per Carton = 300 pcs



Code: TT-09

77" Blood Administration Set Non-Vented

Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 100mm + 70mm (Double Chamber)/
Blood Filter 60mm/Latex Free 'Y' Site/Flow Regulator/Rotating Luer Lock Adaptor with Cap/
2.7mm x 3.9mm PVC Tubing- Length 165 + 15cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 175 pcs



Code: TT-10

64" Blood Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Tube Latex /
Flow Regulator/ Luer Tip Male Adaptor/Needle 18g x 1.5" with Protector/
2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 500 pcs

Option: (TT-10A) Above specification with Bulb Latex instead of Tube Latex

Option: (TT-10B) Above specification with Y Connector & Luer Lock (without Needle)




Code: TT-11

64" Blood Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Tube Latex /
Flow Regulator/ Luer Tip Male Adaptor/Needle 18g x 1.5" with Protector/
2.8mm x 3.8mm PVC Tubing- Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 500 pcs

Option: (TT-11A) Above specification with Y Connector & Luer Lock (without Needle)




Code: TT-13

64" Blood Administration Set Vented

Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Bulb Latex /
Flow Regulator/ Luer Tip Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 500 pcs

Option: (TT-13A) Above specification with Luer Lock instead of Luer Tip




Code: TT-14

64" Blood Administration Set Non-Vented

Non-Vented Spike 20 drops-ml with Guard/Drip Chamber 80mm/Blood Filter /Bulb Latex /
Flow Regulator/ Luer Tip Male Adaptor with Cap/2.8mm x 3.8mm PVC Tubing-
Length 150cm/EO Sterile/Poly Pack
Quantities per Carton = 450 pcs

Option: (TT-14A) Above specification with Luer Lock instead of Luer Tip



Code: TT-15

101" 'Y' Division Blood Administration Set

2 Nos. Non-Vented Spike 15 drops-ml with Guard-Bevel type/2 Nos. Blood Chamber 100mm/
1 no. Blood Filter 60mm/Latex Free 'Y' Site/PVC 2Way Connector/3 Nos. Flow Regulators/
Rotating Luer Lock Male Adaptor with Cap/3.0mm x 4.1mm PVC Tubing-
Length 2 nos. 50+135+15cm/EO Sterile/Paper Poly Pack
Quantities per Carton = 150 pcs